

TESTING MANUAL

Purpose of this manual:

This manual is designed to guide users through performing the three strabismus detection tests without the use of a connected camera. Instead, users can capture images using their own mobile phone cameras and upload them manually into the appropriate folders for each test.

PROCEDURE:

- Note patient name as **Name_id_dob_date**
- Dob and date will be in this format YYYY-MM-DD if it is 27 june 2005 then it is 2005-06-27
- **Ex:-** Teja Kumar Narayana_12345678_2005-06-27_2025-06-26

***Note:** Please ensure that all captured images are clear and that the patient's full face is clearly visible. Failure to meet these criteria may result in an error.*

1. Preliminary Test

Folders Involved:

- preliminary
- Preliminary_Results

Instructions:

1. Store the captured image for the preliminary test in the **preliminary** folder. Name the image according to the convention specified in the earlier steps.
2. Ensure the image you place is the most **recent** among all images in the **preliminary** folder.
3. After completing the steps above, open your browser and navigate to:
<http://localhost:5000/results>
4. If results are not generated, try the following steps:
 1. Clear all data from both **preliminary** and **Preliminary_Results** folders, or move the existing files to a backup location.
 2. Then, place only the current test image in the **preliminary** folder.
5. If you still do not receive any results, proceed to **Step 4**.

Successful Execution:

If the process runs successfully, the results will be saved in the **Preliminary_Results** folder. The output will include:

- A result image with the same **name** used earlier.
- A corresponding text file with the same **name**.

2. Hirschberg Test

NOTE: The background shouldn't be exposed to a bright setting because Hirschberg test contains use of LED light whose reflection in the iris centre is measured for operating the tests.

Folders Involved:

- HirschbergTestImages
- Hirschberg_Results

Instructions:

1. Store the image captured for the Hirschberg test in the **HirschbergTestImages** folder. Use the same naming convention as mentioned in the previous steps.
2. Ensure the image is the most **recent** among all images in the **HirschbergTestImages** folder.
3. After completing the above steps, open your browser and go to: http://localhost:5000/hirschberg_results
4. If results are not being generated, follow these steps:
 1. Clear all data from both the **HirschbergTestImages** and **Hirschberg_Results** folders, or move the existing files to a backup location.
 2. Place only the current test image in the **HirschbergTestImages** folder.
5. If results still do not appear, proceed to **Step 4**.

Successful Execution:

If the process runs successfully, you will find the results in the **Hirschberg_Results** folder. The output will include:

- A result image with the same **name** used earlier.
- A text file named in the format: **processed_{{patient name}}.txt**

3. 9 Gaze Test:

Folders Involved:

- 9GazeTestImages
- 9GazeResults

gaze_1	gaze_2	gaze_3
gaze_4	gaze_5	gaze_6
gaze_7	gaze_8	gaze_9

Instructions:

1. Capture all 9 images in the exact order listed above.
Do not capture images out of sequence (e.g., gaze_1, then gaze_3). The system will treat the first captured image (regardless of name) as gaze_1. Always capture them sequentially from gaze_1 to gaze_9.
2. Name each image according to the format:
gaze_1.jpg, gaze_2.jpg, ..., gaze_9.jpg
3. Create a folder using the naming convention defined in the earlier steps:
{{FullName_PatientID_DOB}}_{{TestDate}}
Example:
Teja Kumar Narayana_12345678_2005-06-27_2025-06-26
 - **DOB** and test **date** should be in the format: **YYYY-MM-DD**.
4. Place all 9 gaze images into this folder.
5. Move the folder into the **9GazeTestImages** folder.
6. In your browser, go to the following URL:
http://localhost:5000/9gaze_results?name={{FullName_PatientID_DOB}}
Example:
http://localhost:5000/9gaze_results?name=Teja Kumar Narayana_12345678_2005-06-27

Successful Execution:

If successful, the results will be stored in the **9GazeResults** folder under a subfolder named:

{{FullName_PatientID_DOB_TestDate}}_Areal

Example:

Teja Kumar Narayana_12345678_2005-06-27_2025-06-26_Areal